3	
}	
į	
į	

Delphion Intellegibil Mangerty Newwork



IPN Home | Search | Order | Shopping Cart | Login | Site Map | Help





View Images (1 pages) | View INPADOC only

Country:

JP Japan

Kind:

Inventor(s):

TANIDA MOCHIMASA HASEGAWA AKIRA OYA TEIICHI

ANRAKU YASUHIRO

Applicant(s):

TONEN CORP

News, Profiles, Stocks and More about this company

Issued/Filed Dates:

Dec. 12, 1995 / May 30, 1994

Application Number:

JP1994000116731

IPC Class:

C12Q 1/02; C12Q 1/02;

Abstract:

Purpose: To simply and rapidly screen an immunosuppressive substance by preparing two sets of yeast media mutually different in Ca concentrations, adding a specimen sample to one of the sets, inoculating a yeast mutant defective in the activity of tonoplast proton ATPase in the media, and culturing the yeast mutant. Constitution: Two sets of yeast culture media mutually different in the concentrations of calcium are prepared. A specimen sample is added to one of the two sets of the yeast culture media, and the specimen sample is not added to the other. A yeast mutant defective in the activity of tonoplast proton ATPase (e.g.

Saccharomyces cerevisiae DV3TL-B2) is inoculated on the media and subsequently cultured. When a calcium concentration allowing the multiplication of the mutant in the medium containing the specimen sample is lower than a calcium concentration allowing the multiplication of the mutant in the medium not containing the specimen sample, it is estimated that the specimen sample contains an immunosuppressive substance. Thereby, the objective immunosuppressive substance for immunosuppressive agents,

etc., can simply and rapidly be screened.

COPYRIGHT: (C)1995,JPO

Other Abstract Info:

CHEMABS 124(11)143613S CAN124(11)143613S DERABS C96-064852

Other Abstract Info:

DERC96-064852

Foreign References:

(No patents reference this one)

Alternative

Searches

Patent Number Boolean Text

Nominate this invention for the Gallery ...

Browse

Disclosure Bulletin

Privacy | Legal | Gallery | IP Pages | Advertising | FAQ | Contact Us